# Sock Aid


## EN Sock Aid

Please inspect all parts carefully for cracks/breaks prior to use.

### USER INSTRUCTIONS

1. For easier control, tie the ends of the straps together.
2. Gather up to the stocking. Squeeze the sock aid into a "U" shape trough. Pull the sock over the narrow end of the trough without the straps.
3. Slide the gathered sock up into the two notches on the side of the sock aid. The sock opening should not extend beyond the two notches and the sock toe should fit snugly over the narrow end. (see figure 1)
4. The heel of the sock should be underneath the sock aid. Use a little talcum powder at the entrance of the sock aid to help the foot slide into the sock aid more easily.
5. Sit on the edge of the bed or chair. Place the sock aid on the floor. Grasp the straps with both hands and guide the sock aid onto your foot. Pull the straps, guiding the sock aid until the sock is completely on your foot. (see figure 2)
6. Repeat process for other foot.

NB. The black rubber circle is fixed to the sock aid purely to stop the stocking from slipping off when being used.

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

Valentine International Ltd. Taiwan  
8th Floor, No. 149, Sec 2, Ta Tung Road  
221 Hsichih District, New Taipei City  
Taiwan R.O.C

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